# Analysis of CLASH data using a Snakemake pipeline Hyb-CRAC-R


Ignatius Pang
23<sup>rd</sup> April 2021 (Updated)

#### Outline


- Introduction to the CLASH experimental protocol
- Go through the steps of the CLASH pipeline
- Go through some concepts in writing Snakemake pipeline
- How to configure the parameters for a new experiment

#### Cross-linking, ligation and sequencing of hybrids (CLASH)

**CLASH:** method to capture and map RNA-RNA interactions using a RNA-binding protein as bait




#### Overall Workflow



**Fig. 1.** Schematic of CLASH experiment and hyb analysis pipeline.

Demultiplexing and Identification of Unique Reads versus PCR duplicate reads



### What is Snakemake (in my own words)?

- Workflow management toolkits (e.g. Nextflow, ruffus, Snakemake)
- Automation like a recipe to run the pipeline, with restart checkpoints

- Define a list of Input and Output files (nodes in the network)
- Define files using "Wildcards"
- Create a network of dependencies
- Run jobs based on available files and dependencies

#### Example of Wildcards

- Suits = ["Diamond", "Club", "Heart", "Spade"]
- Ranks = ["A", "2", "3", "4", "5", "6", "7", "8", "9", "10", "J", "Q", "K"]

- All\_cards= expand("\{rank\}\_of\_\{suit\}.card"), suit= Suits, rank= Ranks),
- {suit} and {rank} are the wildcards
- What if we add {colour}? There are no black diamond?
  - Merge rank and color (e.g. {suit\_colour}) or
  - Use python lookup tables to avoid wrong combinations

### **Example Snakemake Command**

```
Number of parallel jobs that can run at the same time
snakemake -j 12 \
--printshellcmds \
                           Print to STDOUT
                                Your snakemake file
--snakefile jw bc1.smk \
--cluster-config cluster.json \
                                     Json file defines the Katana PBS parameters
--latency-wait 20 --use-envmodules \
                                                20 secs wait between jobs, and use PBS module load system
--cluster "qsub -N {cluster.N} -l
                                          The Katana PBS 'qsub' command, with parameters from JSON file
nodes={cluster.nodes}:ppn={cluster.ppn},walltime={cluster.walltime},mem={clu
ster.mem} -M {cluster.email} -m ae -j oe"
```


We will focus on how we can adjust the file jw\_bc1.smk to perform the analysis.

#### run\_bbmerge run\_dmplx run\_bbduk run\_pcr\_de\_dup run\_hyb\_detect run\_hyb\_merge run\_pyCRAC run\_sgr2pws run annotation get\_max\_reads get\_dataset\_stats calc\_hyb\_p\_val run\_fastqc run\_sep\_unique\_multimapped run\_multiqc

## Rule Graph

```
snakemake \
   --snakefile jw_bc1.smk \
   --rulegraph \
   | dot -Tpdf \
   > rulegraph.pdf
```

## Directed Acyclic Graph (DAG)



```
* A directed acyclic graph of all the steps
snakemake \
--snakefile jw_bc1.smk \
--dag \
| dot -Tpdf \
> dag.pdf
```

## Parameters we need to adjust - Input FASTA Files

## Input Samples Parameter

```
OUTPUT_DIR="/srv/scratch/treelab/igy/Data/CLASH_JKD6009/JW_BC1/Processed_Data"
FILES_R1=["JKD60091-rncHTF-1-2-3-ATCACG_S1_L001_R1_001.fastq.gz",
"JKD60091-rncHTF-1-2-3-ATCACG S1 L002 R1 001.fastq.gz",
"JKD60091-rncHTF-1-2-3-ATCACG_S1_L003_R1_001.fastq.gz",
"JKD60091-rncHTF-1-2-3-ATCACG S1 L004 R1 001.fastq.gz" ]
FILES_R2=["JKD60091-rncHTF-1-2-3-ATCACG_S1_L001_R2 001.fastq.gz",
"JKD60091-rncHTF-1-2-3-ATCACG S1 L002 R2 001.fastq.gz",
"JKD60091-rncHTF-1-2-3-ATCACG_S1_L003_R2_001.fastq.gz",
                                                                             Updated
"JKD60091-rncHTF-1-2-3-ATCACG S1 L004 R2 001.fastq.gz"]
```

FASTQ\_DIR="/srv/scratch/treelab/JJT\_fastq/staph\_CLASH\_may\_2017/W0N3452\_20170510\_Rnc\_Pr

## Parameters we need to adjust - Input Barcode File

```
NNNTAAGC L5Aa_BC1_control_1
NNNATTAGC L5Ab_BC1_rnc-HTF_1 JW_BC1_CLASH_barcodes.tab
NNNGCGCAGC L5Ac_BC1_rnc-HTF-2
NNNCACTAGC L5Bc_BC1_rnc-HTF_3
```

#### Other Parameters We Need To Adjust

```
## Genomic parameters
GENOME DIR="/srv/scratch/treelab/genome/saus/JKD6009"
GTF FILE DIR = os.path.join(GENOME DIR, "Hyb pyCRAC")
pyCRAC_GTF_FILE = os.path.join(GTF_FILE_DIR, "GCF_900607245.1_JKD6009_genomic_merge_features_simple_annot_checked_edited.gtf")
GENOME SIZE = os.path.join(GENOME DIR, "Sa JKD6009.genome")
NOVO_INDEX = os.path.join(GENOME_DIR, "Sa_JKD6009.nix")
HYB_DB_NAME = "Sa_JKD6009"
RNA TYPE PRIORITIES = os.path.join(GTF FILE DIR, "rna type priority.tab")
## Source Directories
SOURCE DIR= os.path.join( BASE DIR, "Source" )
R_SCRIPTS_DIR=os.path.join(SOURCE_DIR, "Demultiplex", "scripts" )
ANNOTATION SCRIPT DIR = os.path.join(SOURCE DIR, "Jai Script", "Annotation")
STATS_SCRIPT_DIR = os.path.join(SOURCE_DIR, "Jai_Script", "Hyb_stats" )
PYCRAC SRC DIR = os.path.join(SOURCE DIR, "Jai Script", "pyCRAC")
## Location of CRAC installation environment
CRAC_ENVIRONMENT_ACTIVATE="/srv/scratch/z3371724/my_python_dir/my_python_env/bin/activate"
```

The pyCRAC environment needs to replaced with a new one that you've installed yourself. (see file StaphCLASH2020/Source/Jai\_Script/using\_pyCRAC\_on\_katana.sh)

## Replace the CRAC\_ENVIRONMENT\_ACTIVATE parameter containing the code which activate python virtualenv environment for pyCRAC scripts:

./srv/scratch/z3371724/my\_python\_dir/my\_python\_env/bin/activate

Updated

# Within the script CRAC\_pipeline\_SE.py, replace this line in the 'cmd' parameter

```
def trimAndCollapse(inputfile,outputfile):
    """ Removes the last random nucleotide from the forward read, collapses the data and then splits it again into fasta files """
    cmd = '/srv/scratch/z3371724/StaphCLASH2020/Source/Jai_Script/pyCRAC/TrimNucs.py    -n 1 --addtoheader -f '%s' | pyFastqDuplicateRemover.py -o '%s''
    logger.info(cmd)
    os.system(cmd)
    #os.system("touch %s" % outputfile)
```

- Replace the 'cmd' parameter within the "trimAndCollapse" function of the script CRAC\_pipeline\_SE.py
- Use the absolute path of the TrimNucs.py script within the GitHub repository (to override and not use the original pyCRAC script)
- For reason not yet determined, the script does not seem to work unless we use an absolute path.
- Would be good to fix this so we do not need to use the absolute path.



Begin With the End in Mind

- The "rule\_all" is like a giant checklist for Snakemake
- Only make files that are not already exist in the directories
- Many files that have a pattern are defined using the wildcards. For example, the {barcode} wildcard.

#### List Out All the Relevant Files in Rule All (part 1)

```
rule all:
    input:
        sample_r1 = expand(os.path.join( FASTQ_DIR, "{fastq}"), fastq=FILES_R1),
        sample_r2 = expand(os.path.join( FASTQ_DIR, "{fastq}"), fastq=FILES_R2),
        linker = LINKER_FILE,
        all_fa = os.path.join(OUTPUT_DIR, "Merged", "all.fq"),
        ihist = os.path.join(OUTPUT_DIR, "Merged", "ihist.txt"),
        barcodes_file = BARCODE_FILE,
        fastq_file= expand( os.path.join(OUTPUT_DIR, "Demultiplexed", "all_{barcode}.fastq"), barcode=BARCODE_NAMES),
        trimmed_files = expand( os.path.join(OUTPUT_DIR, "Trimmed", "all_{barcode}.trimmed.pyCRAC.fq"), barcode=BARCODE_NAMES),
```

There are other rules here but not shown here as it is too long.

Updated

### List Out All the Relevant Files in Rule All (part 2)

Input for BBMerge

```
sample_r1 = expand(os.path.join( FASTQ_DIR, "{fastq}"), fastq=FILES_R1),
sample_r2 = expand(os.path.join( FASTQ_DIR, "{fastq}"), fastq=FILES_R2),
```

Output from BBMerge and Input for Demultiplexing

```
all_fa = os.path.join(OUTPUT_DIR, "Merged", "all.fq"),
```

Output from Demultiplexing and Input for Trimming

```
fastq_file= expand( os.path.join(OUTPUT_DIR, "Demultiplexed", "all_{barcode}.fastq"), barcode=BARCODE_NAMES),
```

Output from Trimming and Input for PCR deduplication

Updated

```
trimmed_files = expand( os.path.join(OUTPUT_DIR, "Trimmed",
    "all {barcode}.trimmed.pyCRAC.fq"), barcode=BARCODE NAMES),
```

There are other rules in 'rule all' but not shown here, it would be too long.

#### run bbmerge run\_dmplx run bbduk run\_pcr\_de\_dup run hyb detect run hyb merge run pyCRAC run\_sgr2pws run annotation get max reads get dataset stats calc\_hyb\_p\_val run fastqc run multiac run\_sep\_unique\_multimapped

### 1. Merge Pairedend Reads

- BBMerge
- Provide adapter sequence in a FASTA file
- Keep the Merged and Unmerged reads and concatenate them into a single file (Not throw away data)

# Defining Input and Output Files Within the Bbmerge Rule

```
rule run_bbmerge:
  input:
    sample_r1 = expand(os.path.join( FASTQ_DIR, "{fastq}"), fastq=FILES_R1),
    sample_r2 = expand(os.path.join( FASTQ_DIR, "{fastq}"), fastq=FILES_R2),
    linker = LINKER FILE
  output:
    all_fa = os.path.join(OUTPUT_DIR, "Merged", "all.fq"),
    ihist = os.path.join(OUTPUT DIR, "Merged", "ihist.txt"),
  shell:
    \mathbf{H} \cdot \mathbf{H} \cdot \mathbf{H}
   cat {input.sample_r1} > $TMPDIR/sample_r1.gz
   cat {input.sample_r2} > $TMPDIR/sample_r2.gz
```

And then the rest of the codes is the run bbmerg using \$TMPDIR/sample\_r1.gz and \$TMPDIR/sample\_r2.gz

#### run bbmerge run\_dmplx run bbduk run\_pcr\_de\_dup run hyb detect run hyb merge run pyCRAC run annotation run\_sgr2pws get max reads get dataset stats calc\_hyb\_p\_val run fastqc run multigc run sep unique multimapped


## 2. Demultiplexing

- Use pyCRAC's pyBarcodeFilter.py
- Notes to install pyCRAC on Katana using python virtuaenv available
- Barcode file

NNNTAAGC L5Aa\_BC1\_control\_1
NNNATTAGC L5Ab\_BC1\_rnc-HTF\_1
NNNGCGCAGC L5Ac\_BC1\_rnc-HTF-2
NNNCACTAGC L5Bc\_BC1\_rnc-HTF\_3

# Defining Input and Output Files Within the Demultiplexing Rule

```
rule run_dmplx:
    input:
    all_fa = os.path.join(OUTPUT_DIR, "Merged", "all.fq"),
    barcodes_file = BARCODE_FILE,
    output:
    os.path.join(OUTPUT_DIR, "Demultiplexed", "all_others.fastq"),
    expand( os.path.join(OUTPUT_DIR, "Demultiplexed", "all_{barcode}.fastq"), barcode=BARCODE_NAMES),
    os.path.join(OUTPUT_DIR, "Demultiplexed", "all_barcode_statistics.txt"),
    os.path.join(OUTPUT_DIR, "Demultiplexed", "all_random_nucleotide_statistics.txt"),
```



## 3. Reads Trimming & QA

#### **Trimming**

- Trim away poor quality sequence regions or adapters
- Note additional parameters in the Snakemake file

#### QA

- Run FASTQC
- Run MultiQC to view all fastQC output at the same time
- Check for over-represented sequences (could be adapters not cleaned properly)

#### run bbmerge run\_dmplx run bbduk run\_pcr\_de\_dup run hyb detect run hyb merge run pyCRAC run\_sgr2pws run annotation get max reads get dataset stats calc\_hyb\_p\_val run fastqc run multiac run sep unique multimapped

### 4. PCR De-duplication

- The hybrid reads were tagged with barcodes prior to PCR amplification
- Each barcode contain two sections

>L5Ab\_BC1\_rnc-HTF\_1 NNNATTAGC

- Random triplet nt (helps identify and remove PCR duplicates)
- ACAATTAGC
- GTCATTAGC

PCR duplicates of the

ACAATTAGC

same read

 5 – 7 nt of known barcode sequence (tell apart different samples)

#### run bbmerge run\_dmplx run bbduk run\_pcr\_de\_dup run hyb detect run hyb merge run pyCRAC run\_sgr2pws run annotation get max reads get dataset stats calc hyb p val run fastqc run multiac run sep unique multimapped

## 5 – 7. Analyze Hybrids

- Use the Hyb pipeline to identify the two-halves of the hybridized read
- Blast of each of the two halves of the hybrid sequence
- Find genomic location of each halve
- Multiple reads representing the same RNA-RNA interactions were merged
- Assign each halve to transcriptome annotations (sRNA, 5'UTR, 3'UTR, coding region)

#### Annotation Step Requires GTF file

- Require one GTF file containing all annotations for all RNA class
- RNA\_TYPE\_PRIORTIES parameter points to the file "rna\_type\_priority.tab". This file list in order of priority of annotation to help resolve overlapping features (see table on the right)
- Directory location of the GTF file is perthe GTF\_FILE\_DIR parameter: /srv/scratch/treelab/genome/saus/JKD6009 /Hyb\_pyCRAC

rank	rna_type
1	CDS
2	5UTR
3	3UTR
4	rRNA
5	tRNA
6	CRISPR
7	pseudogene
8	riboswitch
9	Rnase_P_RNA
10	RNA_termometer
11	SRP_RNA
12	tmRNA
13	operon

Updated

#### run bbmerge run\_dmplx run bbduk run\_pcr\_de\_dup run\_hyb\_detect run hyb merge run\_pyCRAC run sgr2pws run annotation get max reads get dataset stats calc hyb p val run fastqc run multiac run sep unique multimapped

## 8-11. Get Max Reads Depth

- Use the pyCRAC pipeline to get the maximum read depth at each of the two genomic location of the interaction RNAs (run\_pyCRAC, run\_sgr2pws, get\_max\_reads)
- Also creates a file called "datasets.txt" that gives us the total number of reads for each barcode (get\_dataset\_stats)

#### run bbmerge run\_dmplx run bbduk run\_pcr\_de\_dup run hyb detect run hyb merge run pyCRAC run\_sgr2pws run annotation get max reads get dataset stats calc hyb p val run fastqc run multiac run sep unique multimapped

#### 12-13. Calculate p-value

- Uses a binomial distribution to calculate the p-value for each RNA-RNA interaction
- g\_x, g\_y are the maximum read depths of the two RNA
- p\_x\_y=2\*g\_x\*g\_y
- pdf\_x\_y = p\_x\_y/total\_reads
- Hyb\_count = Number of hybrid reads
- dbinom(x=hyb\_count, size= total\_reads, prob = pdf\_x\_y)
- Separate hybrids that are uniquely mapped to the genome with those that are mapped to multiple locations in the genome

# Merge p-values of RNA-RNA interaction from biological replicate experiments

- merge\_p\_values.Rmd
- Uses Fisher's method to combine the p-value (based on a Chi-square distribution)
- Also perform Benjamini-Hochberg FDR correction
- Run independently outside of Snakemake pipeline
- Merge results from multiple Snakemake runs for different biological replicate experiments
- Adjust within R script the following parameters:
- list\_of\_files <- list(</li>
   file.path( hyb\_p\_val\_dir, "DM\_MiSeq", "all\_NNNCACTAGC\_L5Bc\_rnc-HTF\_1.unique.hyb.p\_values.tab" ), file.path(x) etc....)
- names(list\_of\_files ) <- c( "DM\_MiSeq\_Bc", "DM\_NextSeq\_Ba" etc...)</li>

#### How to debug

- Use the --dry-run flag:
  - Simulate all the steps without running them, this checks all the dependencies
  - snakemake --dry-run -j 12 --snakefile jw\_bc1.smk
- Use the --printshellcmds flag:
  - Print all the STDOUT from running each of the rules, useful for debugging
  - Combine with PBS –j oe, -m ae flags to save the output in log file (e.g. job\_name.o12321)
  - Snakemake --printshellcmds -j 12 --snakefile jw\_bc1.smk
- Job killed in the middle of action, and is locked out from a restart
  - snakemake -j 12 --unlock --printshellcmds --snakefile jw\_bc1.smk
- Consider printing out the Rule Graph or the Directed Acyclic Graph

### A quick dry-run demo (successful run)

snakemake --dry-run --printshellcmds -j 12 --snakefile jw\_bc1\_test\_dry\_run.smk

Dry-run completed successfully without fail

```
treelab/igy/Data/CLASH JKD6009/JW BC1/Test Dry Run/Hyb Annot/all NNNTAAGC L5Aa BC1 control 1.multimapped.hyb.p values.tab, /srv/scratch/treelab/ig
 /Data/CLASH JKD6009/JW BC1/Test Dry Run/Hyb Annot/all NNNATTAGC L5Ab BC1 rnc-HTF 1.multimapped.hyb.p values.tab, /srv/scratch/treelab/igy/Data/CLA
   jobid: 0
lob counts:
               jobs
       count
               all
               calc hyb p val
               get dataset stats
               get max reads
               run annotation
               run bbduk
               run bbmerge
               run dmplx
               run fastqc
               run hyb detect
               run hyb merge
               run multiqc
               run pcr de dup
               run pyCRAC
               run sep unique multimapped
               run sgr2pws
This was a dry-run (flag -n). The order of jobs does not reflect the order of execution.
[z3371724@katana1 Script Per Dataset]$
```

## A quick dry-run demo (problem with code)

- I removed the following line from the run\_bbmerge, output declarations:
  - all\_fa = os.path.join(OUTPUT\_DIR, "Merged", "all.fq"),

The run failed with error messages produced.

```
[z3371724@katana1 Script_Per_Dataset]$ snakemake --dry-run --printshellcmds -j 12 --snakefile jw_bc1_test_dry_run.smk

Building DAG of jobs...

MissingInputException in line 59 of /srv/scratch/treelab/igy/Projects/StaphCLASH2020/Source/Demultiplex/Script_Per_Dataset/common_rules_dry_run.smk:

Missing input files for rule run_dmplx:

/srv/scratch/treelab/igy/Data/CLASH_JKD6009/JW_BC1/Test_Dry_Run/Merged/all.fq
```

#### GitHub Resources

https://github.com/lgnatiusPang/Hyb-CRAC-R

https://github.com/lgnatiusPang/TermPick

https://github.com/lgnatiusPang/Snakemake Demo

### A Tip for Writing Snakemake Pipeline

- "Look on every exit as being an entrance somewhere else."
- Tom Stoppard, Rosencrantz and Guildenstern Are Dead

"Look on every output file as being an input file to something else"

## Acknowledgements



- Stuart McKellar
- Sander Granneman

Tim Stinear's Lab

Tree Lab